Device of the Month: Multi-Track Angiographic Catheter

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
Examples of Multi-catheter approaches

- for the determination of pulmonary resistance, the simultaneous reading of pulmonary and aortic pressures is useful. An MT catheter can be positioned in the pulmonary artery and then pushed off the wire. The wire can then be positioned in the aorta (figure 7) and there an MT catheter can be pushed off. The end-hole catheter used for the manipulation of the guidewire can then be used in order to obtain blood samples and pressures and atrial and venous level. This more complex multi-catheter approach is an example for setups in which pharmacological testing is planned. Any variation of catheter positions and combinations can be imagined. (FIGURE 7)

- The guidewire is first positioned and an MT catheter is advanced over the wire and then kicked off the end of the wire. The wire is then positioned elsewhere and another MT catheter is advanced on the newly positioned wire. (icon)

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
Figure 5 The guidewire is positioned in the left pulmonary branch in a patient with tricuspid atresia with transposition of the great arteries who underwent banding of the pulmonary artery. Despite the complicated route of the wire (inferior vena cava — right atrium — atrial septal defect — left atrium — left ventricle — transposed pulmonary artery) the Multi-Track Angio catheter could be advanced easily into the pulmonary artery (A). The pressure gradient over the banding could be measured by a controlled pullback procedure and angiography could be performed just proximal to the band without catheter recoil (B).

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