



Multi-Snare®



INTRADE Interventional Devices

pfm

Application: Loop catching

Application: Lateral catching

Applications

- ▶ Positioning of grafts and implants.
- ▶ Grasping of attached catheters
- ▶ Retrieval of catheters, misplaced implants or other foreign bodies
- ▶ During cross-over procedures

Product features

- ▶ Multi-Snare's dual-plane system increases grasping efficiency offering loop and lateral catching possibility
- ▶ Variable loop diameter for precise sizing of loop
- ▶ Radiopaque loop, platinum coiled
- ▶ Flexible nitinol core wire provides high tensile strength
- ▶ Low profile introducer
- ▶ Low friction with 1:1 torque, during application

Order Information



Multi-Snare® and Multi-Snare® Micro are packaged as a set including:

- ▶ Snare
- ▶ Snare introducer sheath (separate pouch)
- ▶ Snare introduction aid
- ▶ Torquer

Multi-Snare® Sets

REF	Snare		Introducer	
	Diameter	Length	Size	Length
147305	5 mm (5-8)	125 cm	4 F	105 cm
147310	10 mm (10-15)	125 cm	4 F	105 cm
147315	15 mm (15-20)	125 cm	5 F	105 cm
147320	20 mm (20-30)	125 cm	6 F	105 cm
147330	30 mm (30-40)	125 cm	6 F	105 cm

Multi-Snare® single on request

Multi-Snare® Micro Sets

REF	Snare		Introducer	
	Diameter	Length	Size	Length
147302	2 mm (2-3)	175 cm	3 F	150 cm
147304	4 mm (4-6)	175 cm	3 F	150 cm

Multi-Snare® single on request

Introducer Sheath (long)

REF	Sheath Size	Inner Diameter	Required Guide Wire	Introducer Length
147350	3 F	0.021" (0,6 mm)	0.021"	150 cm
147415	4 F	0.035" (1,0 mm)	0.035"	105 cm
147515	5 F	0.044" (1,2 mm)	0.038"	105 cm
147615	6 F	0.055" (1,4 mm)	0.038"	105 cm

Radiopaque tip marker
Low friction catheter with shapeable tip

Certified by DIN EN ISO 9001/EN ISO 13485
pfm - Produkte für die Medizin AG
 Wankelstraße 60 · D-50996 Köln, Germany
 Phone +49 (0) 2236 9641-0 · Fax +49 (0) 2236 9641-20
 customerservice@pfm-ag.com · www.pfm-ag.com



Exclusive Distributor for Europe:

genius
A pfm Company